

Escalade® Acetabular System



## Escalade® Acetabular Shell

- True hemisphere design
- Titanium plasma spray coating
- Designed to maximize head size while maintaining bearing thickness
- Combination taper and snap locking mechanism
- Anti-rotational tabs are flush with the shell on all liners
- +0 liner face truly flush with the shell

## Escalade® Polyethylene

- GUR 1020 resin utilized for superior mechanical properties
- 75kGy gamma irradiated crosslinking process
- Proprietary thermal treatment to eliminate free radicals
- EtO sterilized to maintain crosslinked integrity

## Femoral Head Compatibility

- Ortho Development CoCr heads: 26mm, 28mm, 32mm, 36mm and 40mm
- Biolox® Delta Ceramic heads: 28mm, 32mm, 36mm and 40mm

### Escalade® Compatibility Table

		SHELL DIAMETER														
		42	44	46	48	50	52	54	56	58	60	62	64	66	68	70
HEAD DIAMETER	26	26/42 +3	26/44 +1													
	28		28/44 +3	28/46 +1	28/48	28/50	28/52	28/54	28/56	28/58	28/60	28/62	28/64-70			
	32				32/48 +3	32/50 +1	32/52	32/54	32/56	32/58	32/60	32/62	32/64-70			
	36						36/52 +3	36/54 +1	36/56	36/58	36/60	36/62	36/64-70			
	40								40/56 +3	40/58 +1	40/60	40/62	40/64-70			
<b>HOODED LINERS AVAILABLE IN THESE SIZES</b>																



Ortho Development Corporation  
odev.com

12187 S Business Park Drive  
Draper, UT 84020  
801-553-9991/fax 801-553-9993