

Pivot Bipolar® and Solitude® Unipolar





Pivot Bipolar®

Simple one-step assembly interface between the femoral head and the compression molded liner requires no instrumentation and reduces O.R. time

- 75° range of motion
(12/14 taper and a 28mm femoral head)
- Available in sizes 38mm-60mm in 1mm increments
- Cobalt Chrome heads in 22mm and 28mm



Solitude® Unipolar

Simple intraoperative assembly

- Available in sizes 38mm-60mm in 1mm increments
- Trunion adapters available in -3, +0, +3, +6 and +9 offset options