


Vusion® OS

The Vusion® OS system is designed for versatility and procedural efficiency in posterior lumbar fusions. The system offers intuitive inserter options to accommodate different preferences and techniques while providing a secure implant-inserter interface. The implant's bullet nose allows for distraction during insertion and its convex geometry matches the disc space. In addition, Vusion® OS implants offer:

- Large graft windows to allow for maximum graft placement
- Posteriorly angled teeth designed to provide a secure fit without migration or subsidence
- Tantalum markers for easy intraoperative identification



FOOTPRINT

		9x20mm		9x25mm		9x30mm		11x20mm		11x25mm		11x30mm		11x35mm	
		Parallel	6°	Parallel	6°	Parallel	6°	Parallel	6°	Parallel	6°	Parallel	6°	Parallel	6°
HEIGHT	7mm														
	8mm														
	9mm														
	10mm														
	11mm														
	12mm														
	13mm														
	14mm														
15mm															
16mm															